


Chondroitin Sulfate Methacryloyl (ChSMA)

Product component

Item	character	Package Size	Notes
A: ChSMA	White spongy	0.5 g/bottle	Keep in dark
B: Photoinitiator LAP	White powder	0.025 g/bottle	

This instruction applies to EFL-ChSMA-001



ChSMA molecular structure

Product introduction

Chondroitin sulfate methacryloyl (ChSMA) is a double bond modified chondroitin sulfate and can be quickly photo-crosslinked and cured into gel through UV and visible light in the presence of a photoinitiator. Due to the portable cross-linking method and good biocompatibility, ChSMA-based material systems have been widely used in many biomedical research fields, including: osteoarthritis treatment, joint cartilage repair, skull repair, etc. ChS is rich in carboxyl groups and hydroxyl groups that are easily modified and can be used to build a variety of biomaterials, such as nanodrug carriers and bioadhesives for tumor diagnosis and treatment.

Applications

3D Cell culture, biological 3D printing, tissue engineering, etc.

Storage

Dry kit: room temperature, 3 months; 4°C, 12 months; -20°C, 18 months. **Sterile solution:** 4°C (in dark), 7 days; -20°C (in dark), 6 months. **Please note that repeated freezing and thawing of the solution will affect the performance of the product, so it is best to prepare it when using it.**



企业微信公众号
扫描右侧二维码
获取更多信息

period of validity

The date of manufacture is shown in the package.

Solution preparation

1. Prepare 0.25% (w/v) standard solution of initiator

(1) Add 10ml PBS into the brown bottle containing initiator LAP (containing 0.025g LAP);

(2) Heat and dissolve the solution in a water bath at 40-50°C for 15 minutes, shaking several times.

The LAP standard solution can be stored for 12 months at 4°C in dark.

2. Prepare ChSMA solution (4-10% (w/v) is recommended)

(1) Take the required mass of ChSMA into the centrifugal tube/glass bottle/beaker;

(2) Add the initiator standard solution into the above container;

(3) Dissolve the solution in room temperature for 30 minutes, protected in dark, shaking several times;

(4) Sterilize the ChSMA solution with a 0.22μm sterile needle filter, keep in dark.

Suggestions for 2D cell culture

➤ Inject ChSMA solution into the well plate;

(96-well plate: 50-100 μL/ well, 48-well plate: 100-300 μL/ well, 24-well plate: 300-500 μL/ well)

➤ Irradiate the wells with 405 nm light for 10-30 seconds to gelate, the gel strength can be adjusted by the time and intensity of the light;

➤ Add medium to the wells to cover the gel. Place the well plate in a 37°C incubator for 5 minutes. And then wash the sample and remove the medium;

➤ Add the cell suspension to the well plate. Change medium, observe, and photograph according to experimental design. (No special requirements for operation procedures).

Suggestions for 3D cell culture

➤ Cells were collected and resuspended in ChSMA solution to prepare the cell suspension;

➤ Add cell suspension into the well plates;

(96-well plate: 50-100 μL/ well, 48-well plate: 100-300 μL/ well, 24-well plate:



企业微信公众号
扫描右侧二维码
获取更多信息

300-500 μ L/ well)

- Irradiate the wells with 405 nm light for 10-30 seconds to gelate, the gel strength can be adjusted by the time and intensity of the light;
- Add medium to the wells. Place the plate in a 37°C incubator for 5 minutes. And then wash the sample and remove the medium;
- Add fresh medium and incubate for a long time. Change medium, observe, and photograph according to experimental design. (No special requirements for operation procedures).

Tips: Do not look directly at the light source.



企业微信公众号
扫描右侧二维码
获取更多信息