

Porous Gelatin Methacryloyl

Porous GelMA

Product component

Item	character	Package Size	Notes
Porous GelMA hydrogel	White spongy	0.6 g/ bottle x2	Keep in dark

This instruction applies to EFL-GM-PR-001/002

Storage

Dry kit: room temperature, 3 months; 4°C, 12 months; -20°C, 18 months. **Sterile solution:** 4°C (in dark), 7 days; -20°C (in dark), 6 months. **Please note that repeated freezing and thawing of the solution will affect the performance of the product, so it is best to prepare it when using it.**

period of validity

The date of manufacture is shown in the package.

Required materials

- EFL-GM-PR series porous GelMA hydrogel products^{EFL}
- EFL-LS-1601 series 405nm curing light source equipment^{EFL}
- PBS (1X)
- Constant temperature magnetic stirring water baths
- 0.22 μm Sterile needle filters
- 10-50 mL Sterile centrifuge tubes
- 10 mL Syringes
- 1-5 mL Pipettes & tips



企业微信公众号
扫描右侧二维码
获取更多信息

Operation steps (3D cell culture)

Steps	Title	Materials	Processes
1	Prepare solution	<ul style="list-style-type: none"> ➤ EFL-GM-PR series porous GelMA hydrogel products ➤ PBS (1X) ➤ Pipette guns ➤ constant temperature magnetic stirring water baths 	<ol style="list-style-type: none"> 1) Add appropriate amount of PBS to the porous GelMA bottle; Recommended concentration for porous GelMA is 6-8% w/v. (bottle contains magnetic rotor and 0.6g porous GelMA product). 2) Preparation of porous GelMA hydrogel precursor solution by magnetic stirring at 37°C in a water bath keeping in dark for 1h (important step).
2	Sterilise solution	<ul style="list-style-type: none"> ➤ 0.22 μm Sterile needle filters ➤ Constant temperature water baths ➤ Sterile centrifuge tubes 	<ol style="list-style-type: none"> 1) Sterilise the above solution immediately with a sterile 0.22μm needle filter (to prevent gelation when the temperature drops) and store at 37°C in dark. Note: If it is not possible to use it all at once, store in the refrigerator for a short period of time (< 7 days). Before the next use, redissolve at 37°C and shake for 20-30 seconds to homogenise the material.
3	Mixe cells		<ol style="list-style-type: none"> 1) Collecte cells. 2) Resuspension of cells with sterile precursor solution (multiple blowing or shaking).



企业微信公众号
扫描右侧二维码
获取更多信息

4	Cure GelMA	<ul style="list-style-type: none"> Pipettes & tips EFL-LS-1601 series 405nm curing light source equipment 	<ol style="list-style-type: none"> Add precursor solution to the well plate; 96-well plate: 50-100 μL/ well, 48-well plate: 100-300 μL/ well, 24-well plate: 300-500 μL/ well. Leave at room temperature for 2 min. Irradiation with EFL-LS-1601 series 405nm light source to cure hydrogels, Irradiation times are shown in the table below: <table border="1"> <thead> <tr> <th>Models</th><th>6%</th><th>7%</th><th>8%</th></tr> </thead> <tbody> <tr> <td>001</td><td rowspan="2">16~18s</td><td rowspan="2">13~16s</td><td rowspan="2">8~10s</td></tr> <tr> <td>002</td></tr> </tbody> </table> 	Models	6%	7%	8%	001	16~18s	13~16s	8~10s	002
Models	6%	7%	8%									
001	16~18s	13~16s	8~10s									
002												
5	Wash samples	<ul style="list-style-type: none"> Pipettes & tips 	<ol style="list-style-type: none"> Add medium and incubate at 37°C for 5 minutes. Remove the medium. 									
6	Culture cell		<ol style="list-style-type: none"> Change medium, observe and photograph according to experimental design. 									

Operation steps (2D cell culture)

The main steps of 2D cell culture are the same as those of 3D culture, except that step 3 mixing cells. Once the sample has been washed in step 5, the cells are ready to be grown on the surface.



企业微信公众号
扫描右侧二维码
获取更多信息